

Product Datasheet

CD105 antibody (FITC) (orb44186)

Description

Mouse Monoclonal to CD105.

Reactivity

Human, Porcine

Conjugation

FITC

Tested

FC

Applications
Immunogen

Recombinant Vaccinia virus containing the human CD105 (L-isoform) cDNA.

Target

CD105

Preservatives

Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide

Storage

Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.

Note

For research use only

Application notes

Flow cytometry: The reagent is designed for analysis of human blood cells using 20 µl reagent / 100 µl of whole blood or 106 cells in a suspension. The content of a vial (2 ml) is sufficient for 100 tests.

Isotype

Mouse IgG2a

Clonality

Monoclonal

Purity

Purified antibody is conjugated with fluorescein isothiocyanate (FITC) under optimum conditions and unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.

Uniprot ID

P17813

Entrez


2022

Dilution Range

Flow cytometry: The reagent is designed for analysis of human blood cells using 20 µl reagent / 100 µl of whole blood or 106 cells in a suspension. The content of a vial (2 ml) is sufficient for 100 tests.

Expiration Date

12 months from date of receipt.



Separation of HUVEC cells (red-filled) f...