

Product Datasheet

Human PZP ELISA Kit (orb441863)

Description	Human Pregnancy Zone Protein (PZP) Sandwich ELISA Kit
Reactivity	Human
Range	62.5-4000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P20742
Sensitivity	32.7pg/mL
Entrez	5858
Expiration Date	Please enquire.