

Product Datasheet

Mouse QSOX1 ELISA Kit (orb441867)

Description	Mouse Quiescin Q6 Sulfhydryl Oxidase 1 (QSOX1) Sandwich ELISA Kit
Reactivity	Mouse
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q8BND5
Sensitivity	34pg/mL
Entrez	104009
Expiration Date	Please enquire.