

Product Datasheet

Human RAB5A ELISA Kit (orb441870)

Description	Human RAB5A, Member RAS Oncogene Family (RAB5A) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P20339
Sensitivity	0.122ng/mL
Entrez	5868
Expiration Date	Please enquire.