

Product Datasheet

Human Rac1 ELISA Kit (orb441872)

Description	Human Ras Related C3 Botulinum Toxin Substrate 1 (Rac1) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P63000
Sensitivity	0.060ng/mL
Entrez	5879
Expiration Date	Please enquire.