

Product Datasheet

Rat RANKL ELISA Kit (orb441886)

Description	Rat Receptor Activator Of Nuclear Factor Kappa B Ligand (RANKL) Sandwich ELISA Kit
Reactivity	Rat
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q9ESE2
Sensitivity	5.3pg/mL
Entrez	117516
Expiration Date	Please enquire.