

Product Datasheet

Rat RANK ELISA Kit (orb441887)

Description	Rat Receptor Activator Of Nuclear Factor Kappa B (RANK) Sandwich ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	F1M8Z6
Sensitivity	0.055ng/mL
Entrez	498206
Expiration Date	Please enquire.