

Product Datasheet

Human RARg ELISA Kit (orb441898)

Description	Human Retinoic Acid Receptor Gamma (RARg) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P13631
Sensitivity	0.064ng/mL
Entrez	5916
Expiration Date	Please enquire.