

Product Datasheet

Bovine RETN ELISA Kit (orb441937)

Description	Bovine Resistin (RETN) Sandwich ELISA Kit
Reactivity	Bovine
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q762I5
Sensitivity	34pg/mL
Entrez	369020
Expiration Date	Please enquire.