

Product Datasheet

Rat RGS6 ELISA Kit (orb441946)

Description	Rat Regulator Of G Protein Signaling 6 (RGS6) Sandwich ELISA Kit
Reactivity	Rat
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P49801
Sensitivity	0.118ng/mL
Entrez	54295
Expiration Date	Please enquire.