

Product Datasheet

Porcine RLN ELISA Kit (orb441962)

Description	Porcine Relaxin (RLN) competitive inhibition ELISA Kit
Reactivity	Porcine
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P01348
Sensitivity	8.36pg/mL
Entrez	396891
Expiration Date	Please enquire.