

Product Datasheet

Rat RNASE3 ELISA Kit (orb441968)

Description	Rat Ribonuclease A3 (RNASE3) Sandwich ELISA Kit
Reactivity	Rat
Range	62.5-4000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P70709
Sensitivity	28.73pg/mL
Entrez	192264
Expiration Date	Please enquire.