

Product Datasheet

Bovine ANG ELISA Kit (orb441970)

Description	Bovine Angiogenin (ANG) Sandwich ELISA Kit
Reactivity	Bovine
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P10152
Sensitivity	6.1pg/mL
Entrez	777597
Expiration Date	Please enquire.