

Product Datasheet

Rat RNLS ELISA Kit (orb441985)

Description	Rat Renalase (RNLS) Sandwich ELISA Kit
Reactivity	Rat
Range	6.25-400ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q5U2W9
Sensitivity	2.45ng/mL
Entrez	361751
Expiration Date	Please enquire.