

Product Datasheet

Human RORg ELISA Kit (orb441990)

Description	Human RAR Related Orphan Receptor Gamma (RORg) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P51449
Sensitivity	0.071ng/mL
Entrez	6097
Expiration Date	Please enquire.