

Product Datasheet

Human RPL13A ELISA Kit (orb441993)

Description	Human Ribosomal Protein L13A (RPL13A) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P40429
Sensitivity	0.058ng/mL
Entrez	23521
Expiration Date	Please enquire.