

Product Datasheet

Human RPS6Kb1 ELISA Kit (orb441998)

Description	Human Ribosomal Protein S6 Kinase Beta 1 (RPS6Kb1) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P23443
Sensitivity	0.058ng/mL
Entrez	6198
Expiration Date	Please enquire.