

Product Datasheet

Mouse S100A10 ELISA Kit (orb442018)

Description	Mouse S100 Calcium Binding Protein A10 (S100A10) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P08207
Sensitivity	0.113ng/mL
Entrez	20194
Expiration Date	Please enquire.