

Product Datasheet

Rat S100A4 ELISA Kit (orb442029)

Description	Rat S100A4 ELISA Kit
Reactivity	Rat
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P05942
Sensitivity	34pg/mL
Entrez	24615
Expiration Date	Please enquire.