

Product Datasheet

Human S100A7 ELISA Kit (orb442034)

Description	Human S100A7 ELISA Kit
Reactivity	Human
Range	125-8000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P31151
Sensitivity	51pg/mL
Entrez	6278
Expiration Date	Please enquire.