

Product Datasheet

Mouse S1PR1 ELISA Kit (orb442053)

Description	Mouse Sphingosine 1 Phosphate Receptor 1 (S1PR1) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	008530
Sensitivity	0.115ng/mL
Entrez	13609
Expiration Date	Please enquire.