

Product Datasheet

Mouse SCARA3 ELISA Kit (orb442069)

Description	Mouse Scavenger Receptor Class A Member 3 (SCARA3) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.781-50ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q8C850
Sensitivity	0.33ng/mL
Entrez	219151
Expiration Date	Please enquire.