

## Product Datasheet

### Human SCARA5 ELISA Kit (orb442070)

|                          |                                                                                             |
|--------------------------|---------------------------------------------------------------------------------------------|
| <b>Description</b>       | Human Scavenger Receptor Class A<br>Member 5 (SCARA5) Sandwich ELISA Kit                    |
| <b>Reactivity</b>        | Human                                                                                       |
| <b>Range</b>             | 0.156-10ng/mL                                                                               |
| <b>Form/Appearance</b>   | Traditional                                                                                 |
| <b>Storage</b>           | All reagents should be kept according to the<br>instructions on the packaging upon receipt. |
| <b>Note</b>              | For research use only                                                                       |
| <b>Application notes</b> | tissue homogenates and other biological<br>fluids.                                          |
| <b>Sample Types</b>      | tissue homogenates and other biological<br>fluids.                                          |
| <b>Uniprot ID</b>        | <b>Q6ZMJ2</b>                                                                               |
| <b>Sensitivity</b>       | 0.057ng/mL                                                                                  |
| <b>Entrez</b>            | <b>286133</b>                                                                               |
| <b>Expiration Date</b>   | Please enquire.                                                                             |