

Product Datasheet

Human SCARB1 ELISA Kit (orb442071)

Description	Human Scavenger Receptor Class B Member 1 (SCARB1) Sandwich ELISA Kit
Reactivity	Human
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q8WTV0
Sensitivity	14.3pg/mL
Entrez	949
Expiration Date	Please enquire.