

Product Datasheet

Human SCCA1 ELISA Kit (orb442072)

Description	Human Squamous Cell Carcinoma Antigen 1 (SCCA1) Sandwich ELISA Kit
Reactivity	Human
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P29508
Sensitivity	34pg/mL
Entrez	6317
Expiration Date	Please enquire.