

Product Datasheet

Rat SCGN ELISA Kit (orb442093)

Description	Rat Secretagogin (SCGN) Sandwich ELISA Kit
Reactivity	Rat
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q6R556
Sensitivity	24.7pg/mL
Entrez	306942
Expiration Date	Please enquire.