

Product Datasheet

Rat SCNN1a ELISA Kit (orb442095)

Description	Rat Amiloride Sensitive Sodium Channel Subunit Alpha (SCNN1a) Sandwich ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P37089
Sensitivity	0.059ng/mL
Entrez	25122
Expiration Date	Please enquire.