

Product Datasheet

Rat SDC4 ELISA Kit (orb442108)

Description	Rat Syndecan 4 (SDC4) Sandwich ELISA Kit
Reactivity	Rat
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P34901
Sensitivity	0.133ng/mL
Entrez	24771
Expiration Date	Please enquire.