

Product Datasheet

Human SEMA3B ELISA Kit (orb442141)

Description	Human Semaphorin 3B (SEMA3B) Sandwich ELISA Kit
Reactivity	Human
Range	1.56-100ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q13214
Sensitivity	0.59ng/mL
Entrez	7869
Expiration Date	Please enquire.