

Product Datasheet

Human SEMA3E ELISA Kit (orb442143)

Description	Human Semaphorin 3E (SEMA3E) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	O15041
Sensitivity	0.116ng/mL
Entrez	9723
Expiration Date	Please enquire.