

Product Datasheet

Human SEMA3F ELISA Kit (orb442145)

Description	Human Semaphorin 3F (SEMA3F) Sandwich ELISA Kit
Reactivity	Human
Range	46.875-3000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q13275
Sensitivity	19.25pg/mL
Entrez	6405
Expiration Date	Please enquire.