

Product Datasheet

Human SEMA4B ELISA Kit (orb442147)

Description	Human Semaphorin 4B (SEMA4B) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q9NPR2
Sensitivity	0.071ng/mL
Entrez	10509
Expiration Date	Please enquire.