

Product Datasheet

Human SEMA5A ELISA Kit (orb442149)

Description	Human Semaphorin 5A (SEMA5A) Sandwich ELISA Kit
Reactivity	Human
Range	0.781-50ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q13591
Sensitivity	0.27ng/mL
Entrez	9037
Expiration Date	Please enquire.