

Product Datasheet

Human SEMA7A ELISA Kit (orb442150)

Description	Human Semaphorin 7A (SEMA7A) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	075326
Sensitivity	0.059ng/mL
Entrez	8482
Expiration Date	Please enquire.