

Product Datasheet

Human SFRP5 ELISA Kit (orb442176)

Description	Human Secreted Frizzled Related Protein 5 (SFRP5) Sandwich ELISA Kit
Reactivity	Human
Range	1.56-100ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q5T4F7
Sensitivity	0.60ng/mL
Entrez	6425
Expiration Date	Please enquire.