

Product Datasheet

Mouse SFRP5 ELISA Kit (orb442177)

Description	Mouse Secreted Frizzled Related Protein 5 (SFRP5) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.25-16ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q9WU66
Sensitivity	0.09ng/mL
Entrez	54612
Expiration Date	Please enquire.