

Product Datasheet

Rat SGLT2 ELISA Kit (orb442184)

Description	Rat Sodium/Glucose Cotransporter 2 (SGLT2) Sandwich ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P53792
Sensitivity	0.051ng/mL
Entrez	64522
Expiration Date	Please enquire.