

Product Datasheet

Human SHBG ELISA Kit (orb442192)

Description	Human Sex Hormone Binding Globulin (SHBG) Sandwich ELISA Kit
Reactivity	Human
Range	62.5-4000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P04278
Sensitivity	27.8pg/mL
Entrez	6462
Expiration Date	Please enquire.