

Product Datasheet

Filamin A antibody (orb4422)

Description	Filamin A antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human FLNA around the phosphorylation site of Ser2151 (AP(p-S)VA)
Target	phospho-Filamin A (Ser2151)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	291 kDa
Uniprot ID	P21333
Expiration Date	12 months from date of receipt.