

Product Datasheet

Human SIRPg ELISA Kit (orb442213)

Description	Human Signal Regulatory Protein Gamma (SIRPg) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q9P1W8
Sensitivity	0.057ng/mL
Entrez	55423
Expiration Date	Please enquire.