

Product Datasheet

Mouse SIRT1 ELISA Kit (orb442215)

Description	Mouse Sirtuin 1 (SIRT1) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q923E4
Sensitivity	0.055ng/mL
Entrez	93759
Expiration Date	Please enquire.