

Product Datasheet

Rat SIRT3 ELISA Kit (orb442220)

Description	Rat Sirtuin 3 (SIRT3) Sandwich ELISA Kit
Reactivity	Rat
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	B2RZ31
Sensitivity	0.088ng/mL
Entrez	293615
Expiration Date	Please enquire.