

Product Datasheet

Mouse SLC30A7 ELISA Kit (orb442239)

Description	Mouse Solute Carrier Family 30 Member 7 (SLC30A7) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q9JKN1
Sensitivity	0.053ng/mL
Entrez	66500
Expiration Date	Please enquire.