

Product Datasheet

Rat SNCa ELISA Kit (orb442280)

Description	Rat Synuclein Alpha (SNCa) Sandwich ELISA Kit
Reactivity	Rat
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P37377
Sensitivity	5.4pg/mL
Entrez	29219
Expiration Date	Please enquire.