

Product Datasheet

Rat SOCS1 ELISA Kit (orb442286)

Description	Rat Suppressors Of Cytokine Signaling 1 (SOCS1) Sandwich ELISA Kit
Reactivity	Rat
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q9QX78
Sensitivity	0.112ng/mL
Entrez	252971
Expiration Date	Please enquire.