

Product Datasheet

Human SOCS3 ELISA Kit (orb442288)

Description	Human Suppressors Of Cytokine Signaling 3 (SOCS3) Sandwich ELISA Kit
Reactivity	Human
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	O14543
Sensitivity	6.4pg/mL
Entrez	9021
Expiration Date	Please enquire.