

Product Datasheet

Human SOST ELISA Kit (orb442307)

Description	Human Sclerostin (SOST) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q9BQB4
Sensitivity	0.122ng/mL
Entrez	50964
Expiration Date	Please enquire.