

## Product Datasheet

### Canine SPA ELISA Kit (orb442314)

|                          |                                                                                          |
|--------------------------|------------------------------------------------------------------------------------------|
| <b>Description</b>       | Canine Surfactant Associated Protein A (SPA)<br>Sandwich ELISA Kit                       |
| <b>Reactivity</b>        | Canine                                                                                   |
| <b>Range</b>             | 1.875-120ng/mL                                                                           |
| <b>Form/Appearance</b>   | Traditional                                                                              |
| <b>Storage</b>           | All reagents should be kept according to the instructions on the packaging upon receipt. |
| <b>Note</b>              | For research use only                                                                    |
| <b>Application notes</b> | serum, plasma, lung lavage fluid, tissue homogenates and other biological fluids.        |
| <b>Sample Types</b>      | serum, plasma, lung lavage fluid, tissue homogenates and other biological fluids.        |
| <b>Uniprot ID</b>        | <b>Q6RXL1</b>                                                                            |
| <b>Sensitivity</b>       | 0.62ng/mL                                                                                |
| <b>Entrez</b>            | <b>407213</b>                                                                            |
| <b>Expiration Date</b>   | Please enquire.                                                                          |