

Product Datasheet

Porcine SPA ELISA Kit (orb442317)

Description	Porcine Surfactant Associated Protein A (SPA) Sandwich ELISA Kit
Reactivity	Porcine
Range	3.125-200ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, lung lavage fluid and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, lung lavage fluid and other biological fluids.
Uniprot ID	P49874
Sensitivity	1.16ng/mL
Entrez	397503
Expiration Date	Please enquire.