

Product Datasheet

Rat SPA ELISA Kit (orb442318)

Description	Rat Surfactant Associated Protein A (SPA) Sandwich ELISA Kit
Reactivity	Rat
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, lung lavage fluid, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, lung lavage fluid, tissue homogenates and other biological fluids.
Uniprot ID	P08427
Sensitivity	6.9pg/mL
Entrez	24773
Expiration Date	Please enquire.